



Email : sales@laboratoryinstrumentindia.com

Product Name :
Jackets with Hood Snaps

Product Code :
LBNY-00017-11000808



Description :

Material: Polyester, PVC
Color: Yellow
No. of Pockets: 2
Size (English) Chest: 66 in.
Size (Metric) Chest: 167.6 cm
Length (English) Sleeve: 36.5 in.
Length (Metric) Sleeve: 92.7 cm.
Size: 3X-Large

Technical Specification :

Jackets with Hood Snaps



Laboratory Instrument India