



Email : sales@laboratoryinstrumentindia.com

Product Name :
CTI Spunbonded Polypropylene Shoe Covers

Product Code :
LBNY-00017-11000812



Description :

Made from spun bond polypropylene
Lightweight material offers comfort and breathability, while maintain snug, full coverage
Individually inspected for stitching, sizing, materials and workmanship
Lot traceability
Strong and durable shoe covers are made from a lightweight material that offers breathability and comfort.
Size: X-Large
Color: Blue.

Technical Specification :

CTI Spunbonded Polypropylene Shoe Covers



Laboratory Instrument India