



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
CTI Spunbonded Polypropylene Shoe Covers

**Product Code :**  
LBNY-00017-11000812



**Description :**

Made from spun bond polypropylene  
Lightweight material offers comfort and breathability, while maintain snug, full coverage  
Individually inspected for stitching, sizing, materials and workmanship  
Lot traceability  
Strong and durable shoe covers are made from a lightweight material that offers breathability and comfort.  
Size: X-Large  
Color: Blue.

**Technical Specification :**

CTI Spunbonded Polypropylene Shoe Covers



**Laboratory Instrument India**