



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Disposable Coveralls with Hood and Boot

**Product Code :**  
LBNY-00017-11000814



**Description :**

Disposable Coveralls are designed and tested to protect against blood, blood-borne pathogens, bodily fluids, and many other liquid hazards.

Roomy, seamless shoulder area resists tears and provides more protection in key exposure area

With an attached hood and boots, you're protected from head-to-toe

Blue color makes stains and splashes less conspicuous

Possible applications include emergency response, category 3 black water, tank cleaning, crime scene remediation, and animal research

These coveralls boast elastic wrists and back, and a front zipper with storm flap.

Tear- and abrasion-resistant fabric protects against scrapes and scratches in the toughest of work environments

Fabric helps prevent heat stress with enhanced moisture vapor transmission, so you stay cool and comfortable on the job

Wider, gusseted crotch gives wearer more room to move in all directions without ripping fabric

Elastic back stretches with user for more comfort and better tear resistance.

Fabric tests for protection against blood, bodily fluids and blood borne pathogens.

**Technical Specification :**

Disposable Coveralls with Hood and Boot

---

# Laboratory Instrument India

