



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Die-Cut Aprons

**Product Code :**  
LBNY-00017-11000821



**Description :**

With an economical and lightweight one-piece construction for comfort, these aprons are ideal for critical and controlled environments  
Available in two lengths and three thicknesses.  
Separate ties eliminated to avoid contamination problems

**Technical Specification :**

Die-Cut Aprons



**Laboratory Instrument India**