



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Pocket Zipper Mesh Vest

Product Code :
LBNY-00017-11000822



Description :

Features reinforced front closure for rugged wear and tear.
Polyester mesh construction (100%)
Zipper front closure reinforced with durable webbing for rugged wear/tear
Breathable contrasting trim reduces hot spots
Gusseted radio/utility pocket
Hook and loop front closure (Diameter: 2 in.)
Left chest mic tabs for quick access
Left chest pocket
Reflective piping (Width: 2 in.)
Piping down sides, utility pocket for added low light visibility.
Reflective trim around arms.

Technical Specification :

Single Pocket Zipper Mesh Vest



Laboratory Instrument India

