



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Single Pocket Zipper Mesh Vest

**Product Code :**  
LBNY-00017-11000822



**Description :**

Features reinforced front closure for rugged wear and tear.  
Polyester mesh construction (100%)  
Zipper front closure reinforced with durable webbing for rugged wear/tear  
Breathable contrasting trim reduces hot spots  
Gusseted radio/utility pocket  
Hook and loop front closure (Diameter: 2 in.)  
Left chest mic tabs for quick access  
Left chest pocket  
Reflective piping (Width: 2 in.)  
Piping down sides, utility pocket for added low light visibility.  
Reflective trim around arms.

**Technical Specification :**

Single Pocket Zipper Mesh Vest



**Laboratory Instrument India**

