



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
ESD, Anti-Static Counter Lab Jackets

**Product Code :**  
LBNY-00017-11000824



**Description :**

Electronic blue counter jacket features a five-gripper front as well as adjustable gripper sleeves with grounding snaps.

Left chest pocket

Two lower pockets.

One-piece collar with notched lapel

Adjustable gripper sleeve with grounding snaps.

**Technical Specification :**

ESD, Anti-Static Counter Lab Jackets



**Laboratory Instrument India**