



Email : sales@laboratoryinstrumentindia.com

Product Name :
ESD, Anti-Static Counter Lab Jackets

Product Code :
LBNY-00017-11000824



Description :

Electronic blue counter jacket features a five-gripper front as well as adjustable gripper sleeves with grounding snaps.

Left chest pocket

Two lower pockets.

One-piece collar with notched lapel

Adjustable gripper sleeve with grounding snaps.

Technical Specification :

ESD, Anti-Static Counter Lab Jackets



Laboratory Instrument India