



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bouffant Dispenser

**Product Code :**  
LBNY-00017-11000826



**Description :**

Made of 1/8 in.-thick (0.3cm) clear acrylic  
Full-length hinge for durability  
Flame treated edges for smoothness  
For bulk dispensing in cleanrooms and laboratories  
Freestanding or wall mountable.  
Extra-large rectangular opening for easy access.

**Technical Specification :**

Bouffant Dispenser



**Laboratory Instrument India**