



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
3-Compartment Dispenser with Front Access

**Product Code :**  
LBNY-00017-11000827



**Description :**

Full-length hinge for durability  
Mounts to wall or free-standing  
1/4 in. clear PETG material - SAFE with IPA  
Sloped lid to maintain HEPA airflow  
Flame treated edges for smoothness.  
Wall bracket included.

**Technical Specification :**

3-Compartment Dispenser with Front Access



**Laboratory Instrument India**