



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automated Digital Pressure Test Kit

Product Code :
LBNY-00017-11000829



Description :

Includes
Fittings, hoses, universal adapter and complete set-up and operating instructions.
Requires
Test is performed automatically, allowing the user to handle other duties while the test is in process
Digital readout for more accurate testing
Large LCD touchscreen allows easy test criteria selection and data input
Air supply with a 1/4 inch MNPT fitting.
Improves the efficiency and accuracy of gas-tight garment tests for both ASTM and CE methods
Test data can be copied to the provided USB stick via the built-in USB port for record keeping.

Technical Specification :

Automated Digital Pressure Test Kit



Laboratory Instrument India