



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automated Digital Pressure Test Kit

**Product Code :**  
LBNY-00017-11000829



**Description :**

Includes

Fittings, hoses, universal adapter and complete set-up and operating instructions.

Requires

Test is performed automatically, allowing the user to handle other duties while the test is in process

Digital readout for more accurate testing

Large LCD touchscreen allows easy test criteria selection and data input

Air supply with a 1/4 inch MNPT fitting.

Improves the efficiency and accuracy of gas-tight garment tests for both ASTM and CE methods

Test data can be copied to the provided USB stick via the built-in USB port for record keeping.

**Technical Specification :**

Automated Digital Pressure Test Kit



**Laboratory Instrument India**