



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens Navy Industrial Pant

**Product Code :**  
LBNY-00017-11000830



#### **Description :**

The strong and durable twill fabric has superior soil-release and color-retention capabilities so you won't be spending any time working out stains or figuring out what wash cycle is right, either.

Primary Closure: Zipper

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

These easy-care properties mean our Industrial Pants are always ready to go to work.

Waistband: Inner lined for body and shape, folder set band with outlet

Technology provides breathable comfort, lasting colors and stain resistance with an even softer feel

Black Fabric resists fading wash after wash

Fabric: 7.5 oz. Twill

Blend: 65% Polyester/35% Cotton

Finish: Wrinkle Resistant

Closure: Heavy-duty brass ratcheting zipper, button closure

Pocket: Two slack-style front pockets, two set-in hip pockets, left has button closure, darts over hip pockets for better fit

Fit Details: Traditional work pant fit; straight leg opening

Gender: Male

Materials: Twill, Cotton/Polyester

Requires minimum ironing thanks to a wrinkle-resistant finish

Synthetic blend pocketing and waistband trim

Silhouette: Mens Style - traditional easy fit.

#### **Technical Specification :**

---

Mens Navy Industrial Pant



**Laboratory Instrument India**