



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Chill-Its Evaporative Cooling Vest

**Product Code :**  
LBNY-00017-11000831



**Description :**

Evaporative Cooling Vest provides cooling relief and comfort once soaked in cold water.  
Water-repellent liner and black cotton-poly elastic trim  
Lightweight, durable, and washable  
Soak in water 2-5 minutes to activate  
V-neck with zipper closure  
Quilted nylon outer layer  
Activated, polymer-embedded fabric inside  
Remains hydrated for up to 4 hours  
Mesh side panels offer improved ventilation and stretch.

**Technical Specification :**

Chill-Its Evaporative Cooling Vest



**Laboratory Instrument India**