



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chill-Its Evaporative Cooling Vest

Product Code :
LBNY-00017-11000831



Description :

Evaporative Cooling Vest provides cooling relief and comfort once soaked in cold water.
Water-repellent liner and black cotton-poly elastic trim
Lightweight, durable, and washable
Soak in water 2-5 minutes to activate
V-neck with zipper closure
Quilted nylon outer layer
Activated, polymer-embedded fabric inside
Remains hydrated for up to 4 hours
Mesh side panels offer improved ventilation and stretch.

Technical Specification :

Chill-Its Evaporative Cooling Vest



Laboratory Instrument India