



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Unisex Short-Length Lab Jackets

**Product Code :**  
LBNY-00017-11000833



**Description :**

Short-length lab coats provide AAMI PB70 Level 2 fluid resistance against incidental splash of blood, body fluids, aerosols, and all aqueous liquids. Suitable for cleanroom use.

Primary Applications:

Clinical Laboratory

Phlebotomy

Anatomic Pathology

Clinical Pathology

Specimen Collection

360-degree protection front, shoulders, and back

Fluid resistant fabric, AAMI PB70 Level 2 Barrier

Blood Collection, Testing

High barrier protection against strike through from incidental splash of blood

No cotton fabric used on shoulders or back as it would not meet OSHA 29 if splashed with blood or OPIM traveling rearward either side of the employee's head

Overall air permeability of coat front and back is 18 CFM

Stain fabric and cuffs (blood, bleach and ALL aqueous reagents and liquids)

Anti-static fabric prevents false instrument readings, damage to computers

Low lint fabric prevents sample contamination, ideal for molecular testing

Durable construction, 3 pockets, double stitched

Short-length 31.5 in.

Provides all height wearers protection to ankles when worn in combination with style Unisex Drawstring Pant

High air permeability (35 CFM) is used to make the back of coat for wearer comfort

Stainless snap closures

---

Durable cuffs for over gloving, bleach safe all colors  
Hang loop outside back of collar  
Side vents for pant pocket access.

**Technical Specification :**

Unisex Short-Length Lab Jackets



**Laboratory Instrument India**