



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PF Zero Unisex Lab Coats

**Product Code :**  
LBNY-00017-11000834



**Description :**

Breathable lab coats provide Level 2 fluid resistance and feature heavy-duty stainless-steel snap closures and hip side vents.

Primary Applications:

Clinical Laboratory

Phlebotomy

Specimen Collection

Anatomic Pathology

Clinical Pathology

Blood Collection, Testing

Dental

Human and Animal Research

Provides 360-degree front and back, plus shoulder protection

Overall Air Permeability 35 CFM

Stainless snap closures

Cuffs for over gloving

Bleach safe all colors

Hang loop, outside back of collar

Stain resistant fabric and cuffs (blood, bleach and ALL aqueous reagents and liquids)

Anti-static fabric prevents false instrument readings, damage to computers

Low lint fabric prevents sample contamination, ideal for molecular testing

Durable construction, 3 pockets, double stitched

225 commercial washes tested.

Side vents for pant pocket access

---

Autoclavable fabric.

**Technical Specification :**

PF Zero Unisex Lab Coats



**Laboratory Instrument India**