



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Brilliant Series Police Vest

Product Code :
LBNY-00017-11000835



Description :

Lightweight solid polyester front with Ultra-Cool 100% mesh shoulders and back
Two 1 in. wide high performance reflective material with contrast
Brilliant Trim reflective edging improves low light visibility
Lightweight solid polyester front with 100% mesh shoulders and back
4-point breakaway at shoulders and sides
4-season adjustability allows for accurate sizing in any weather
When not in use, the vest folds into a convenient storage pouch.

Technical Specification :

Premium Brilliant Series Police Vest



Laboratory Instrument India