



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Premium Brilliant Series Police Vest

**Product Code :**  
LBNY-00017-11000835



**Description :**

Lightweight solid polyester front with Ultra-Cool 100% mesh shoulders and back  
Two 1 in. wide high performance reflective material with contrast  
Brilliant Trim reflective edging improves low light visibility  
Lightweight solid polyester front with 100% mesh shoulders and back  
4-point breakaway at shoulders and sides  
4-season adjustability allows for accurate sizing in any weather  
When not in use, the vest folds into a convenient storage pouch.

**Technical Specification :**

Premium Brilliant Series Police Vest



**Laboratory Instrument India**