



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Black Series Soft Shell Jacket

Product Code :
LBNY-00017-11000838



Description :

Made with high quality materials to increase durability. Premium Black Series Soft Shell Jacket gives enhanced visibility with reflective trim and piping.

Reflective piping accents add better visibility in low light

Convenient underarm vents for breathability

Waist cinches hidden in front pockets

Durable water-resistant soft-shell material

2 in. wide high performance reflective material

Cuffs feature durable rubber pull tabs for easy adjustments

Black areas offer contrast, style and durability

Pockets:

Left chest media-compatible pocket

Two outside lower slash pockets with zipper

Left arm cell phone pocket.

Technical Specification :

Premium Black Series Soft Shell Jacket

Laboratory Instrument India

