



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ripstop Bomber Jacket

**Product Code :**  
LBNY-00017-11000839



**Description :**

Features waterproof seam sealed construction and fleece body for protection against tough environments.  
Adjustable cuffs with rubber pulls  
Durable ripstop material with waterproof seam sealed construction  
Left and right chest mic tabs  
Tablet pocket with storm flap  
Hood fits into collar or can be removed altogether  
Quilted sleeves for easy arm in/out access  
Reflective material (Width: 2 in.)  
Pockets  
Right chest vertical pocket for wallet  
Left chest radio pocket  
Large vertical tablet (iPad) pocket with zipper inside storm flap  
Tablet can also be accessed without unzipping jacket.  
Left and right outside slash pockets  
Left bicep 2-division pencil/pen pocket with drainage grommets.

**Technical Specification :**

Ripstop Bomber Jacket

---

# Laboratory Instrument India

