



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ripstop Bomber Jacket

Product Code :
LBNY-00017-11000839



Description :

Features waterproof seam sealed construction and fleece body for protection against tough environments.
Adjustable cuffs with rubber pulls
Durable ripstop material with waterproof seam sealed construction
Left and right chest mic tabs
Tablet pocket with storm flap
Hood fits into collar or can be removed altogether
Quilted sleeves for easy arm in/out access
Reflective material (Width: 2 in.)
Pockets
Right chest vertical pocket for wallet
Left chest radio pocket
Large vertical tablet (iPad) pocket with zipper inside storm flap
Tablet can also be accessed without unzipping jacket.
Left and right outside slash pockets
Left bicep 2-division pencil/pen pocket with drainage grommets.

Technical Specification :

Ripstop Bomber Jacket

Laboratory Instrument India

