



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Laminated Polypropylene Water-Repellant Shoe
Cover with Non-Skid AQ Sole

Product Code :
LBNY-00017-11000842



Description :

Construction: Polypropylene
For Use With (Application): Manufacturing, food service, industrial, medical, clean room, pharmaceutical
Material: Polypropylene
Latex Free: Yes
Size: 2X-Large
Color: White
Sole Type: Non-skid
Dimensions (L x W): 15.25 x 5.12 in. (38.73 x 13 cm).

Technical Specification :

Safety Laminated Polypropylene Water-Repellant Shoe Cover with Non-Skid AQ Sole



Laboratory Instrument India