



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
High Performance Static Dissipative Cleanroom Frocks

**Product Code :**  
LBNY-00017-11000854



**Description :**

Suitable for exposure to bloodborne pathogens, some hazardous liquids, and any wet general maintenance areas

Features include a snap front, elastic wrists

Lint-free film coated frock for use in critical, ESD, and low particulate clean environments.

Finished fabric contains a proprietary additive making it static dissipative

For use in ISO 6 Cleanroom Class environments

Neatly folded and cleanroom packed

Lint-free film coated frock

Mandarin snap collar, no pockets

Fabric is constructed from a blend of polyethylene and polypropylene

10 pieces per sealed inner polybag, 3 polybags per sealed master bag, 30 pieces per case.

**Technical Specification :**

High Performance Static Dissipative Cleanroom Frocks



**Laboratory Instrument India**

