



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unisex 41 in. (104 cm) Long Open Back Gown- Fluid Resistant Barrier

Product Code :
LBNY-00017-11000857



Description :

Open Back Gown with One Piece Seamless Front and High Neck Coverage provides Wearers Front, Shoulders Level 2 Protection from Incidental splash of Blood, Body Fluids, Aerosols, etc.; Open Back provides Exceptional Comfort to Wearer; Stain Resistant, Anti-static; 225 Commercial Wash Tested.

Wearer Protection Provided

No patch pockets on front for splash of blood, body fluids, aerosols or OPIM to accumulate.

Open Back Design offers protection to shoulders from splash traveling either side of employees face.

One Piece Seamless front protects wearers from neck to below knees, arms, and covers tops of legs with no gaping. Textured Fabric holds fabric into position on top of legs when seated.

Adjustable snaps at neck for complete coverage. Back ties at waist.

Technical Specification :

Unisex 41 in. (104 cm) Long Open Back Gown- Fluid Resistant Barrier



Laboratory Instrument India

