



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Unisex 41 in. (104 cm) Long Closed Back Gown- Fluid Resistant Barrier (AAMI Level 2)

**Product Code :**

LBNY-00017-11000859



**Description :**

Closed back gown provides AAMI Level 2 protection from incidental splash of blood, bodily fluids, and aerosols.

Wearer Protection Provided

Adjustable Snaps at Neck for full Coverage. Back Ties at Waist

No patch pockets on front for splash of blood or OPIM to accumulate.

One Piece Seamless front protects wearers from neck to below knees, arms, and covers tops of legs with no gaping. Textured Fabric holds fabric into position on top of legs when seated.

Closed Back Design offers protection to shoulders and back from splash traveling either side of employee's face.

Primary Applications

Anatomic Pathology • Clinical Pathology • Clinical Laboratory • Phlebotomy • Specimen Collection • Blood Collection, Testing, & Manufacturing • Human & Animal Research • Academia • Dialysis • Dental • Clean Room ISO 5/Class 100 - Medical Device, Pharma, Electronics.

**Technical Specification :**

Unisex 41 in. (104 cm) Long Closed Back Gown- Fluid Resistant Barrier (AAMI Level 2)

---

# Laboratory Instrument India

