



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Level 2 Isolation Gown

**Product Code :**  
LBNY-00017-11000860



**Description :**

Level 2 isolation gown protects from low-risk fluids while providing a secure fit.  
Proprietary over-the-head design with easy tear away removal.  
Elastic wrist with thumb loop to provide a more secure fit under gloves.  
Meets Level 2 requirements for fluid protection.  
Used in low risk situations.

**Technical Specification :**

Level 2 Isolation Gown



**Laboratory Instrument India**