



Email : sales@laboratoryinstrumentindia.com

Product Name :
Level 2 Isolation Gown

Product Code :
LBNY-00017-11000860



Description :

Level 2 isolation gown protects from low-risk fluids while providing a secure fit.
Proprietary over-the-head design with easy tear away removal.
Elastic wrist with thumb loop to provide a more secure fit under gloves.
Meets Level 2 requirements for fluid protection.
Used in low risk situations.

Technical Specification :

Level 2 Isolation Gown



Laboratory Instrument India