



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Dark Navy Industrial Work Pant

Product Code :
LBNY-00017-11000871



Description :

Work pants built of sturdy fabric for hard wear and easy care. Tradesmen swear by the comfort and durability.

Fabric: Body: 8.5 oz. Twill

Blend: Body: 65% Polyester / 35% Cotton

Finish: Wrinkle Resistant, Stain Release

The comprised of heavyweight 8.5 oz. twill.

Number of Pockets: 4

Primary Closure: Hook and Eye

Wash Care: Home Wash

Distinctive tunnel belt loops

Closure: Metal hook and eye

Pocket: Welt back pockets

Waistband: Inner-lined, 2 in.

Product Collection: Flex

Gender: Male.

Technical Specification :

Mens Dark Navy Industrial Work Pant

Laboratory Instrument India

