



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mens Olive Green Industrial Work Pant

**Product Code :**  
LBNY-00017-11000870



**Description :**

Work pants built of sturdy fabric for hard wear and easy care. Tradesmen swear by the comfort and durability. The comprised of heavyweight 8.5 oz. twill.  
Blend: Body: 65% Polyester / 35% Cotton  
Finish: Wrinkle Resistant, Stain Release  
Closure: Metal hook and eye  
Pocket: Welt back pockets  
Number of Pockets: 4  
Primary Closure: Hook and Eye  
Wash Care: Home Wash  
Fabric: Body: 8.5 oz. Twill  
Product Collection: Flex  
Waistband: Inner-lined, 2 in.  
Distinctive tunnel belt loops  
Gender: Male.

**Technical Specification :**

Mens Olive Green Industrial Work Pant

---

# Laboratory Instrument India

