



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Masses Hooked 20gms

**Product Code :**  
LBNY-0003-1010046



**Description :**

Masses Hooked 20gms

**Technical Specification :**

Hook is free to move and hooks on both sides for additional weights to be connected.

Weight is stainless steel silver color with stamped identifier.

This weight has a bored hole through the weight to allow for the moving double hook to sit flat if need be.

Material: Stainless steel.

Weight: 20gms.



**Laboratory Instrument India**