



Email : sales@laboratoryinstrumentindia.com

Product Name :
Men Dark Navy Premium Industrial Cargo Pant

Product Code :
LBNY-00017-11000873



Description :

Cargo pants offering uniform comfort while standing up to industrial laundering. A wrinkle-free polyester/cotton twill with technology resists fading even in dark colors. Reinforced front pockets, metal tack button and wider belt loops add durability.

Fabric: Body: 7.25 oz. Twill
Blend: Body: 65% Polyester / 35% Cotton
Finish: Wrinkle Resistant, Moisture Wicking, Stain Release
Closure: Metal tack button
Number of Pockets: 6
Primary Closure: Button
Wash Care: Industrial Laundry - Heavy Soil
Pocket: Reinforced front pocket bags, Cargo pocket with hidden zipper closure and stitched down cargo pocket flaps, Flap-free welt back pockets.
Fit Details: Sits slightly below the waist, Extra room in seat and thigh, Straight leg
Waistband: Inner-lined, 1.5 in.
Permanent press finish
Gender: Male

Technical Specification :

Men Dark Navy Premium Industrial Cargo Pant



Laboratory Instrument India