



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Men Dark Navy Premium Industrial Cargo Pant

**Product Code :**  
LBNY-00017-11000873



**Description :**

Cargo pants offering uniform comfort while standing up to industrial laundering. A wrinkle-free polyester/cotton twill with technology resists fading even in dark colors. Reinforced front pockets, metal tack button and wider belt loops add durability.

Fabric: Body: 7.25 oz. Twill  
Blend: Body: 65% Polyester / 35% Cotton  
Finish: Wrinkle Resistant, Moisture Wicking, Stain Release  
Closure: Metal tack button  
Number of Pockets: 6  
Primary Closure: Button  
Wash Care: Industrial Laundry - Heavy Soil  
Pocket: Reinforced front pocket bags, Cargo pocket with hidden zipper closure and stitched down cargo pocket flaps, Flap-free welt back pockets.  
Fit Details: Sits slightly below the waist, Extra room in seat and thigh, Straight leg  
Waistband: Inner-lined, 1.5 in.  
Permanent press finish  
Gender: Male

**Technical Specification :**

Men Dark Navy Premium Industrial Cargo Pant



**Laboratory Instrument India**