



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Taffeta Shoe Cover with CSM Soft Sole

Product Code :
LBNY-00017-11000874



Description :

This taffeta shoe cover with CSM Soft Sole is ideal for any cleanroom that does not require ESD. The CSM sole offers flexibility and comfort, while the taffeta fabric will provide a durable particle barrier.

Size/Color Coded Ribbon

Taffeta Fabric

Full Tunneled Elastic Top

Ankle-Length

Soft Tan CSM Sole

Primary Applications: Cleanroom Pharmaceutical, Biotechnological, Medical Device Manufacturing.

Technical Specification :

Protective Wear Taffeta Shoe Cover with CSM Soft Sole



Laboratory Instrument India