



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sterile Cleanroom Sleeves

**Product Code :**  
LBNY-00017-11000876



**Description :**

Sterile Cleanroom Sleeves are made of a comfortable fabric that protects against various liquids and particulates.

Each sleeve is individually packaged and sterilized to a 10<sup>-6</sup> sterility assurance level (SAL) preventing workplace contamination and cross-contamination.

Sterile sleeves are made of a comfortable fabric that protects against various liquids and particulates.

This also contains raw edges that can create and release additional particulates. Sleeves are tear-resistant and made in an ISO-certified facility to ensure the highest quality.

The low-linting fabric helps reduce the release of unwanted particulates to prevent product damage and workplace contamination. Sleeves are designed with tunnelized elastic, providing a more secure closure and stronger seams.

Generous 21.5 in. length provides excellent arm coverage

Sterile-Don packaging protects against contamination during donning

Gamma Guard helps keep the wearer cooler vs many competitive fabrics

The tear-resistant fabric ensures your safety in challenging environments.

**Technical Specification :**

Sterile Cleanroom Sleeves

---

# Laboratory Instrument India

