



Email : sales@laboratoryinstrumentindia.com

Product Name :
Light Blue SMS Lab Coat with 3 Pockets, Knit Wrists

Product Code :
LBNY-00017-11000877



Description :

SMS lab coat protects against dust, dirt, grime, and dry particulates in nonhazardous settings. Fabric also helps protect against light liquid splashes and sprays. Fabric is breathable to help keep you cool, so you can work comfortably over long hours. With three roomy pockets, you can store eyeglasses, pens, and other small instruments to keep with you on the go. 1 left chest pocket and 2 hip pockets
Possible applications include pharmaceuticals, labs, healthcare, controlled environments, food processing, and more.
Light blue SMS lab coat protects against dust, dirt, grime, and other dry particulates in non-hazardous environments.
Three-layer SMS is lightweight and breathable to help prevent heat stress
Knit wrists and knit collar.

Technical Specification :

Light Blue SMS Lab Coat with 3 Pockets, Knit Wrists



Laboratory Instrument India

