



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Light Blue SMS Lab Coat with 3 Pockets, Knit Wrists

**Product Code :**  
LBNY-00017-11000877



**Description :**

SMS lab coat protects against dust, dirt, grime, and dry particulates in nonhazardous settings. Fabric also helps protect against light liquid splashes and sprays. Fabric is breathable to help keep you cool, so you can work comfortably over long hours. With three roomy pockets, you can store eyeglasses, pens, and other small instruments to keep with you on the go. 1 left chest pocket and 2 hip pockets  
Possible applications include pharmaceuticals, labs, healthcare, controlled environments, food processing, and more.  
Light blue SMS lab coat protects against dust, dirt, grime, and other dry particulates in non-hazardous environments.  
Three-layer SMS is lightweight and breathable to help prevent heat stress  
Knit wrists and knit collar.

**Technical Specification :**

Light Blue SMS Lab Coat with 3 Pockets, Knit Wrists



**Laboratory Instrument India**

