



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
SMS Pleated Bouffant Caps

**Product Code :**  
LBNY-00017-11000878



**Description :**

SMS pleated bouffants are lightweight and breathable  
Durable, three-layer SMS fabric keeps hair securely in place  
Protects hair from contaminants like dust and other particles, while also protecting products from cross-contamination  
Pleated bouffants are lightweight, breathable, and provide an secure particulate barrier.  
Three-layer SMS fabric keeps hair in place, while also keeping worker's cool for extended periods.  
Potential applications include healthcare, pharmaceutical, labs, critical environments, food processing, and more.  
Expandable fit with an elastic closure comfortably accommodates a wide range of hairstyles without pulling when removed.

**Technical Specification :**

SMS Pleated Bouffant Caps



**Laboratory Instrument India**

