



Email : sales@laboratoryinstrumentindia.com

Product Name :
SMS Pleated Bouffant Caps

Product Code :
LBNY-00017-11000878



Description :

SMS pleated bouffants are lightweight and breathable
Durable, three-layer SMS fabric keeps hair securely in place
Protects hair from contaminants like dust and other particles, while also protecting products from cross-contamination
Pleated bouffants are lightweight, breathable, and provide an secure particulate barrier.
Three-layer SMS fabric keeps hair in place, while also keeping worker's cool for extended periods.
Potential applications include healthcare, pharmaceutical, labs, critical environments, food processing, and more.
Expandable fit with an elastic closure comfortably accommodates a wide range of hairstyles without pulling when removed.

Technical Specification :

SMS Pleated Bouffant Caps



Laboratory Instrument India

