



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bouffant Bootie Dispenser

Product Code :
LBNY-00017-11000879



Description :

Dispense booties, bouffants, and other similarly sized apparel with ease.
Versatile design allows for dispensation of a variety of apparel supplies
Top-loading hinged lid for easy re-stocking
The AK Bouffant/Bootie Dispenser features a round, front-based center opening for easy removal of stored apparel and is available in two sizes.
Can be wall-mounted or used free standing
Made of 0.125 in high impact PETG.

Technical Specification :

Bouffant Bootie Dispenser



Laboratory Instrument India