



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bouffant Bootie Dispenser

**Product Code :**  
LBNY-00017-11000879



**Description :**

Dispense booties, bouffants, and other similarly sized apparel with ease.  
Versatile design allows for dispensation of a variety of apparel supplies  
Top-loading hinged lid for easy re-stocking  
The AK Bouffant/Bootie Dispenser features a round, front-based center opening for easy removal of stored apparel and is available in two sizes.  
Can be wall-mounted or used free standing  
Made of 0.125 in high impact PETG.

**Technical Specification :**

Bouffant Bootie Dispenser



**Laboratory Instrument India**