



Email : sales@laboratoryinstrumentindia.com

Product Name :
Guardian Protective Wear Full-Length Raincoats

Product Code :
LBNY-00017-11000881



Description :

Material: Nylon
Closure Type: Snap
Coating: Polyurethane
Size: 4X-Large
Seam Type: Stitched, taped
Color: Yellow
No. of Pockets: 2
Style: Full length
Size (English) Chest: 60 to 62 in.

Technical Specification :

Guardian Protective Wear Full-Length Raincoats



Laboratory Instrument India