



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tablet Oscilloscope (Book Series)

**Product Code :**  
LBNY-0002-620028



**Description :**

Tablet Oscilloscope (Book Series)

**Technical Specification :**

Bandwidth: 200 MHz.  
Rise Time: 1.75ns,  
Sample rate: 2 GS/s(1 CH), Max Wave form Capture  
Rate: 500,000 wfm/s, Input Channel: 2)

Tablet Oscilloscope (Book Series), Tablet Oscilloscope (Book Series) Equipments, Tablet Oscilloscope (Book Series) Tools, Tablet Oscilloscope (Book Series) Tool Kits, Tablet Oscilloscope (Book Series) Manufacturers, Tablet Oscilloscope (Book Series) Suppliers from India, China, Kenya



**Laboratory Instrument India**