



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Guardian Protective Wear Full-Length Raincoats

**Product Code :**  
LBNY-00017-11000882



**Description :**

Material: Nylon  
Closure Type: Snap  
Coating: Polyurethane  
Size: 4X-Large  
Color: Yellow  
No. of Pockets: 2  
Style: Full length  
Seam Type: Stitched, taped  
Size (English) Chest: 60 to 62 in.

**Technical Specification :**

Guardian Protective Wear Full-Length Raincoats



**Laboratory Instrument India**