




Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety SMS Disposable Anti-Static Frocks, No Pockets

Product Code :
LBNY-00017-11000883



Description :

Provide particulate barrier and liquid repellent
Ideal for food processing, laboratories, biohazards, medical, metal surface treating, pharmaceutical
Knit cuffs with elastic thumb loops, six snap closure, knit collar, extra-long knee-length, hanger loop
Knee-length frocks constructed of SMS composite fabric feature knit wrists and elastic thumb loops.
Made from SMS composite fabric
Disposable, without pockets
Highly breathable and comfortable
They provide particulate and liquid barrier suitable for ISO Class 6 cleanrooms.
Low shedding and anti-static to protect sensitive products.

Technical Specification :

Safety SMS Disposable Anti-Static Frocks, No Pockets



Laboratory Instrument India