



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Cool Coveralls with Hood

Product Code :
LBNY-00017-11000884



Description :

Microporous coverall with SMS back creates more comfort through greater breathability in warm work environments.

Great for applications such as Pharmaceutical, Painting, Spray Booth, Industrial manufacturing and Research and clinical facilities.

Breathability in warmer work environments.

Keyguard Microporous Coverall w/ SMS Back

Elastic Wrist, Elastic Ankles

100% Latex Free.

Technical Specification :

Safety Cool Coveralls with Hood



Laboratory Instrument India