



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Masses Slotted

**Product Code :**  
LBNY-0003-1010048



**Description :**

Masses Slotted

**Technical Specification :**

Stainless Steel Slotted Masses Set can be used with hanging balances as well as for calibration and studying pressure, strength, torque in the classroom or lab.

Unit includes 9 weights, and one hanger, each weighing 50g.

Weight: 5 gms.

Material: Stainless Steel



**Laboratory Instrument India**