



Email : sales@laboratoryinstrumentindia.com

Product Name :

Men Lightweight Navy Flame Resistant Mock Turtleneck

Product Code :

LBNY-00017-11000895



Description :

A comfortable and durable CAT 2 AR/FR, mock turtleneck that will keep you warm and protected on the job. Men Lightweight Navy Flame Resistant Mock Turtleneck is a comfortable and durable CAT 2 AR/FR compliant mock turtleneck that will keep you warm and protected on the job.

Hazard: Flash Fire, Arc Flash

Primary Closure: Pull On

Wash Care: Home Wash

Features: Double-layer, rib knit collar sized to protect the neck, Tagless brand label and warning label

Product Collection: Excel FR

Fabric: Flame resistant, 7 oz. (198.4 g/m²) Interlock

Blend: 100% Cotton

Collar: Double-layer, rib knit collar sized to protect the neck

Gender: Unisex

Materials: Cotton

Protection: ATPV Arc Rating: 8.9

Fit Details: Regular Fit.

Technical Specification :

Men Lightweight Navy Flame Resistant Mock Turtleneck

Laboratory Instrument India

