



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Men Navy Flame Resistant Flex Knit Long Sleeve T-shirt

**Product Code :**  
LBNY-00017-11000899



**Description :**

Men Navy Flame Resistant Flex Knit Long Sleeve T-shirt gives you the stretch of your favorite T-shirt while keeping you cool and dry.  
Unique front pocket design with pencil stall  
Ribbed collar  
left sleeve label leaves no question that you are protected on the job site  
Lightweight fabric flexes for added comfort and range of movement  
Fabric: Body: 6.75 oz. (230 g/m<sup>2</sup>;) Flame Resistant Jersey, Contrast Panels: 6.75 oz. (230 g/m<sup>2</sup>)  
Flame Resistant Jersey  
Contrasting pattern on sleeve for an athletic look  
Curved flatlock seams for optimal fit  
Hazard: Flash Fire, Arc Flash  
Wash Care: Home Wash  
Blend: Body: 93% Cotton/7% Spandex Elastane, Contrast Panels: 93% Cotton/7% Spandex Elastane.

**Technical Specification :**

Men Navy Flame Resistant Flex Knit Long Sleeve T-shirt

---

# Laboratory Instrument India

